

PARTIALLY EDENTULOUS CBCT SCAN OF MAXILLA ARCH

1. Patient Positioning:

- i** - Breathe through your nose
- Bite gently on cotton rolls
- Close your eyes and place tongue at the back of the mouth

2. Vital anatomy to be captured:

- i** - Superior: Mid sinus / below orbit
- Inferior: Capture 1-3 mm of mandible dentition



3. CBCT General Settings: Maxl 6cm-8cm, 0.3-0.4 Voxel

4. Export and Upload CBCT Scan

- i** - Export the CBCT scan in a multi-file DICOM format to the desktop
- Right click on the folder, go to "Send To", and select "Compressed (Zipped) Folder"
- Upload the compressed / zip folder to Implant Concierge

5. Ship the diagnostic model to Implant Concierge

- i** - For more information on exporting DICOM and printing a shipping label, please visit: www.implantconcierge.com/support

(866) 977-2228

www.implantconcierge.com

©2017, Implant Concierge, LLC.

All rights reserved. Implant Concierge and Virtual Treatment Plan Coordinator are trademarks of Implant Concierge, LLC.

This product is produced by U.S Patient No. 8.630.869