

---

● **STEP 1**

Patient Positioning:

- Breathe through your nose
- Bite gently on cotton rolls
- Close your eyes and place tongue at the back of the mouth

---

● **STEP 2**

Vital anatomy to be captured:

- Superior: Sinus (below orbit)
- Inferior: Chin and entire mandible



---

● **STEP 3**

CBCT General Settings: Full 13cm-20cm, 0.3-0.4 Voxel

---

● **STEP 4**

Export and Upload CBCT Scan

- Export the CBCT scan in a multi-file DICOM format to the desktop
- Right click on the folder, go to "Send To", and select "Compressed (Zipped) Folder"
- Upload the zip / compressed folder to your Implant Concierge case

---

● **STEP 5**

Ship the diagnostic model to Implant Concierge.

- For more information on exporting DICOM and printing a shipping label, please [www.implantconcierge.com/support](http://www.implantconcierge.com/support)